

Compound T129171

Chemical Properties

CAS No. :

Formula: C₁₆H₂₈O₁₂

Molecular Weight: 412.39

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(2R,3R,4S,5R,6R)-3,5-dihydroxy-2-(hydroxymethyl)-6-(((2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy)oxan-4-yl 2-methylpropanoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4249 mL	12.1244 mL	24.2489 mL
5 mM	0.485 mL	2.4249 mL	4.8498 mL
10 mM	0.2425 mL	1.2124 mL	2.4249 mL
50 mM	0.0485 mL	0.2425 mL	0.485 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481