


## Sahandol

## Chemical Properties

CAS No. :	1452398-07-0
Formula:	C <sub>20</sub> H <sub>24</sub> O <sub>2</sub>
Molecular Weight:	296.4
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sahandol shows antiprotozoal activity against <i>P. falciparum</i> and <i>T. brucei rhodesiense</i> .
Targets(IC50)	Parasite,Antifection

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3738 mL	16.8691 mL	33.7382 mL
5 mM	0.6748 mL	3.3738 mL	6.7476 mL
10 mM	0.3374 mL	1.6869 mL	3.3738 mL
50 mM	0.0675 mL	0.3374 mL	0.6748 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481