


## Trimetozine

## Chemical Properties

CAS No. :	635-41-6
Formula:	C <sub>14</sub> H <sub>19</sub> NO <sub>5</sub>
Molecular Weight:	281.3
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Trimetozine (Sedoxazine) is a sedative that has been marketed in Europe since 1959. It has mild tranquilizing effects and has been used to treat anxiety.
Targets(IC50)	Others,Antibacterial,Antibiotic

## Solubility Information

Solubility	DMSO: 25 mg/mL (88.87 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (7.11 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5549 mL	17.7746 mL	35.5492 mL
5 mM	0.711 mL	3.5549 mL	7.1098 mL
10 mM	0.3555 mL	1.7775 mL	3.5549 mL
50 mM	0.0711 mL	0.3555 mL	0.711 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Shpak VM, Shcheglova Alu. *Pediatrics*. 1968 Aug;47(8):76-7.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481