

Compound T129178

Chemical Properties


CAS No. :

Formula: C₂₅H₃₂O₈

Molecular Weight: 460.52

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

20-hydroxy-4-(hydroxymethyl)-10-methoxy-1,14,19,19-tetramethyl-2,7,18-trioxapentacyclo[11.9.0.0^{3,11}.0^{5,9}.0^{14,20}]docosa-3,5(9),10-triene-8,17-dione is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1715 mL	10.8573 mL	21.7146 mL
5 mM	0.4343 mL	2.1715 mL	4.3429 mL
10 mM	0.2171 mL	1.0857 mL	2.1715 mL
50 mM	0.0434 mL	0.2171 mL	0.4343 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481