


SB-429201

Chemical Properties

CAS No. : 1027971-34-1
 Formula: C₂₈H₂₄N₂O₃
 Molecular Weight: 436.5
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	SB-429201 is an effective, selective inhibitor of HDAC1.
Targets(IC50)	Others,HDAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.291 mL	11.4548 mL	22.9095 mL
5 mM	0.4582 mL	2.291 mL	4.5819 mL
10 mM	0.2291 mL	1.1455 mL	2.291 mL
50 mM	0.0458 mL	0.2291 mL	0.4582 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Hu E, Dul E, Sung CM, Chen Z, Kirkpatrick R, Zhang GF, Johanson K, Liu R, Lago A, Hofmann G, Macarron R, de los Frailes M, Perez P, Krawiec J, Winkler J, Jaye M. Identification of novel isoform-selective inhibitors within class I histone deacetylases. J Pharmacol Exp Ther. 2003 Nov;307(2):720-8. Epub 2003 Sep 15. PubMed PMID: 12975486.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481