


4,6-Di-tert-butyl-3-methylphenol

Chemical Properties

CAS No. :	497-39-2
Formula:	C ₁₅ H ₂₄ O
Molecular Weight:	220.35
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	4,6-Di-tert-butyl-3-methylphenol is a biochemical reagent used in biochemical experiments and pharmaceutical synthesis research.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5382 mL	22.6912 mL	45.3823 mL
5 mM	0.9076 mL	4.5382 mL	9.0765 mL
10 mM	0.4538 mL	2.2691 mL	4.5382 mL
50 mM	0.0908 mL	0.4538 mL	0.9076 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wang, X.; Meng, X.; Shi, L.; Liu, N.; Wang, J. Solid Acid Catalysts in Tert-Butyl Removal Reaction of 4,6-DBMC: Performance Insights and Industrial Implications. Catal. Lett. 2025, 155, 165.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481