


GBR 12935 dihydrochloride

Chemical Properties

CAS No. :	67469-81-2
Formula:	C ₂₈ H ₃₆ Cl ₂ N ₂ O
Molecular Weight:	487.5
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	GBR 12935 dihydrochloride (GBR 12935) can induce the increase of the extracellular levels of dopamine to basal levels.
Targets(IC50)	Dopamine Receptor
Kinase Assay	Recombinant ROCK-I, ROCK-II, PKN, or citron kinase is expressed in HeLa cells as Myc-tagged proteins by transfection using Lipofectamine, and is precipitated from the cell lysates by the use of 9E10 monoclonal anti-Myc antibody coupled to G protein-Sepharose. Recovered immunocomplexes are incubated with various concentrations of [32P]ATP and 10 mg of histone type 2 as substrates in the absence or presence of various concentrations of either Y-27632 or Y-30141 at 30°C for 30 min in a total volume of 30 µL of the kinase buffer containing 50 mM HEPES-NaOH, pH 7.4, 10 mM MgCl ₂ , 5 mM MnCl ₂ , 0.02% Brij 35, and 2 mM dithiothreitol. PKCa is incubated with 5 µM [32P]ATP and 200 µg/mL histone type 2 as substrates in the absence or presence of various concentrations of either Y-27632 or Y-30141 at 30°C for 10 min in a kinase buffer containing 50 mM Tris-HCl, pH 7.5, 0.5 mM CaCl ₂ , 5 mM magnesium acetate, 25 µg/mL phosphatidyl serine, 50 ng/mL 12-O-tetradecanoylphorbol-13-acetate and 0.001% leupeptin in a total volume of 30 µL. Incubation is terminated by the addition of 10 µL of 43 Laemmli sample buffer. After boiling for 5 min, the mixture is subjected to SDS-polyacrylamide gel electrophoresis on a 16% gel. The gel is stained with Coomassie Brilliant Blue, and then dried. The bands corresponding to histone type 2 are excised, and the radioactivity is measured[1].

Solubility Information

Solubility	DMSO: 50 mg/mL (102.56 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (4.1 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0513 mL	10.2564 mL	20.5128 mL
5 mM	0.4103 mL	2.0513 mL	4.1026 mL
10 mM	0.2051 mL	1.0256 mL	2.0513 mL
50 mM	0.041 mL	0.2051 mL	0.4103 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Norlén M, et al. Eur J Pharmacol. 1999 Feb 5;366(2-3):329-32.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481