

Isoginsenoside Rh3

Chemical Properties

CAS No. :	166040-90-0
Formula:	C ₃₆ H ₆₀ O ₇
Molecular Weight:	604.87
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Isoginsenoside Rh3 is a natural product
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6532 mL	8.2662 mL	16.5325 mL
5 mM	0.3306 mL	1.6532 mL	3.3065 mL
10 mM	0.1653 mL	0.8266 mL	1.6532 mL
50 mM	0.0331 mL	0.1653 mL	0.3306 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481