

Resencatinib

Chemical Properties


CAS No. : 2546117-79-5

Formula: C₃₀H₂₉N₇O₃

Molecular Weight: 535.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Resencatinib, a potent tyrosine kinase inhibitor, exhibits antineoplastic activity.
Targets(IC50)	Others,c-Met/HGFR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8671 mL	9.3353 mL	18.6706 mL
5 mM	0.3734 mL	1.8671 mL	3.7341 mL
10 mM	0.1867 mL	0.9335 mL	1.8671 mL
50 mM	0.0373 mL	0.1867 mL	0.3734 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481