


3,4-Dimethoxyphenethylamine

Chemical Properties

CAS No. :	120-20-7
Formula:	C ₁₀ H ₁₅ NO ₂
Molecular Weight:	181.23
Storage:	Store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	3,4-Dimethoxyphenethylamine is a human endogenous metabolite.
Targets(IC50)	Endogenous Metabolite

Solubility Information

Solubility	DMSO: 95 mg/mL (524.2 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.5179 mL	27.5893 mL	55.1785 mL
5 mM	1.1036 mL	5.5179 mL	11.0357 mL
10 mM	0.5518 mL	2.7589 mL	5.5179 mL
50 mM	0.1104 mL	0.5518 mL	1.1036 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Sargent T W , Israelstam D M , Shulgin A T , et al. A note concerning the fate of the 4-methoxyl group in 3, 4-dimethoxyphenethylamine (DMPEA)[J]. Biochemical & Biophysical Research Communications, 1967, 29(1):126-130.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481