


## NH-bis(C1-Boc)

### Chemical Properties

CAS No. :	85916-13-8
Formula:	C <sub>12</sub> H <sub>23</sub> NO <sub>4</sub>
Molecular Weight:	245.32
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



### Biological Description

Description	NH-bis(C1-Boc) is an uncleavable linker used for antibody-drug conjugates (ADC).
Targets(IC50)	ADC Linker, PROTAC Linker

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0763 mL	20.3815 mL	40.7631 mL
5 mM	0.8153 mL	4.0763 mL	8.1526 mL
10 mM	0.4076 mL	2.0382 mL	4.0763 mL
50 mM	0.0815 mL	0.4076 mL	0.8153 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481