

Compound T128196

Chemical Properties


CAS No. :

Formula: C₄₁H₆₆O₁₂

Molecular Weight: 750.97

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2R,4aR,6aS,6bR,10S,12aR,14bS)-10-(((2S,3R,4S,5S)-4,5-dihydroxy-3-(((2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxy)oxan-2-yl]oxy)-2-(hydroxymethyl)-2,6a,6b,9,9,12a-hexamethyl-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14b-icosahydricene-4a-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3316 mL	6.6581 mL	13.3161 mL
5 mM	0.2663 mL	1.3316 mL	2.6632 mL
10 mM	0.1332 mL	0.6658 mL	1.3316 mL
50 mM	0.0266 mL	0.1332 mL	0.2663 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481