


Flavidin

Chemical Properties

CAS No. :	83924-98-5
Formula:	C15H12O3
Molecular Weight:	240.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Flavidin shows very good antioxidant capacity. Flavidin can enhance fluorescent imaging, allowing more sensitive and specific cell labeling in tissues, it should have wide application in molecular detection, providing a general insight into how to optimize simultaneously the behavior of the biomolecule and the chemical probe.
Targets(IC50)	COX,Prostaglandin Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1623 mL	20.8117 mL	41.6233 mL
5 mM	0.8325 mL	4.1623 mL	8.3247 mL
10 mM	0.4162 mL	2.0812 mL	4.1623 mL
50 mM	0.0832 mL	0.4162 mL	0.8325 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481