


Compound PDK0416

Chemical Properties

CAS No. :	95-15-8
Formula:	C ₈ H ₆ S
Molecular Weight:	134.19
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Compound PDK0416, with CAS No. 95-15-8, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Compound PDK0416 provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.4521 mL	37.2606 mL	74.5212 mL
5 mM	1.4904 mL	7.4521 mL	14.9042 mL
10 mM	0.7452 mL	3.7261 mL	7.4521 mL
50 mM	0.149 mL	0.7452 mL	1.4904 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481