


SU11274

Chemical Properties

CAS No. : 658084-23-2
 Formula: C₂₈H₃₀ClN₅O₄S
 Molecular Weight: 568.09
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	SU11274 (Met Kinase Inhibitor)(IC ₅₀ =10 nM) is a specific Met inhibitor. It shows no significant effects on PGDFRβ, EGFR or Tie2.
Targets(IC ₅₀)	Apoptosis,FGFR,CDK,Autophagy,c-Met/HGFR,VEGFR
In vivo	SU11274 inhibits the phosphorylation of key regulators in the PI3K pathway, including AKT, FKHR, and GSK3β. It also inhibits the HGF-dependent phosphorylation of Met and the subsequent HGF-dependent cellular proliferation and activity, with an IC ₅₀ ranging from 1 to 1.5 μM. In BaF3 cells transformed by TPR-MET, the presence of interleukin-3 and treatment with SU11274 suppress cell growth in a dose-dependent manner. However, SU11274 (IC ₅₀ <3 μM) does not inhibit the growth of BaF3 cells transformed by other oncogenic tyrosine kinases, including BCR-ABL, TEL-JAK2, TEL-ABL, and TEL-PDGFR. Treatment with SU11274 also significantly inhibits cell migration in BaF3.TPR-MET cells with inhibitions of 44.8% and 80% at concentrations of 1 μM and 5 μM, respectively. In H69 and H345 cells, which possess functional Met receptors, SU11274 inhibits HGF-induced cell growth with IC ₅₀ values of 3.4 μM and 6.5 μM, respectively. In non-small cell lung cancer (NSCLC) cells expressing c-Met, SU11274 inhibits cell viability with an IC ₅₀ of 0.8-4.4 μM and abrogates the phosphorylation of c-Met and its downstream signaling induced by hepatocyte growth factor.
Kinase Assay	In vitro Met kinase assay: A chimeric protein is constructed containing the cytoplasmic domain of human c-Met fused to Glutathione S-transferase (GST) and expressed in SF9 cells. The c-Met kinase GST-fusion protein is used for an ELISA-based Met biochemical assay using the random copolymer poly(Glu:Tyr) (4:1) immobilized on microtiter plates as a substrate. IC ₅₀ value is determined with various concentrations of SU11274 in a buffer containing 5 μM ATP and 10 mM MnCl ₂ , 50 mM HEPES (pH 7.5), 25 mM NaCl, 0.01% BSA, and 0.1 mM Na orthovanadate. The kinase reaction is performed for 5 minutes at room temperature. The extent of substrate phosphorylation is measured using horseradish peroxidase-conjugated anti-pTyr antibodies.
Cell Research	Cells are exposed to various concentrations of SU11274 in the presence or absence of HGF for 24, 48, and 72 hours. The number of viable cells is determined using the MTT assay or trypan blue exclusion. Cell Cycle and apoptosis are measured by fluorescence-activated cell sorter analysis via propidium iodide staining and Annexin V-positive staining, respectively. (Only for Reference)

Solubility Information

Solubility	H2O: < 1 mg/mL (insoluble or slightly soluble), Ethanol: 2 mg/mL (3.52 mM), Sonication is recommended. DMSO: 48.33 mg/mL (85.07 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 4.83 mg/mL (8.5 mM), Suspension. 10% DMSO+90% Saline: < 4.83 mg/mL (8.5 mM), Lower concentrations may be soluble, but exact solubility limit is unknown. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7603 mL	8.8014 mL	17.6028 mL
5 mM	0.3521 mL	1.7603 mL	3.5206 mL
10 mM	0.176 mL	0.8801 mL	1.7603 mL
50 mM	0.0352 mL	0.176 mL	0.3521 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wang X, et al. Mol Cancer Ther, 2003, 2(11):1085-1092.

Sattler M, et al. Cancer Res, 2003, 63(17), 5462-5469.

Ma PC, et al. Cancer Res, 2005, 65(4), 1479-1488.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481