

## Desthiobiotin-PEG4-propargyl

## Chemical Properties

CAS No. : 1951424-89-7


Formula: C<sub>21</sub>H<sub>37</sub>N<sub>3</sub>O<sub>6</sub>

Molecular Weight: 427.542

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Desthiobiotin-PEG4-propargyl is a PEG-coupled PROteolysis TARgeting Chimeric (PROTAC) linker used in PROTAC synthesis.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands joined by a linker; one ligand targets an E3 ubiquitin ligase, while the other targets the specific protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.339 mL	11.6948 mL	23.3896 mL
5 mM	0.4678 mL	2.339 mL	4.6779 mL
10 mM	0.2339 mL	1.1695 mL	2.339 mL
50 mM	0.0468 mL	0.2339 mL	0.4678 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481