

NHPI-PEG3-C2-NHS ester

Chemical Properties

CAS No. : 2101206-14-6


Formula: C₂₁H₂₄N₂O₁₀

Molecular Weight: 464.42

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	NHPI-PEG3-C2-NHS ester is a noncleavable, three-unit PEG linker used specifically for the synthesis of antibody-drug conjugates (ADCs).
Targets(IC50)	ADC Linker, PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1532 mL	10.7661 mL	21.5322 mL
5 mM	0.4306 mL	2.1532 mL	4.3064 mL
10 mM	0.2153 mL	1.0766 mL	2.1532 mL
50 mM	0.0431 mL	0.2153 mL	0.4306 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481