


Acetildenafil

Chemical Properties

CAS No. :	831217-01-7
Formula:	C ₂₅ H ₃₄ N ₆ O ₃
Molecular Weight:	466.58
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Acetildenafil is an analog of the phosphodiesterase inhibitor sildenafil.
Targets(IC50)	PDE

Solubility Information

Solubility	DMSO: 3.3 mg/mL (7.07 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1433 mL	10.7163 mL	21.4326 mL
5 mM	0.4287 mL	2.1433 mL	4.2865 mL
10 mM	0.2143 mL	1.0716 mL	2.1433 mL
50 mM	0.0429 mL	0.2143 mL	0.4287 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Venhuis B J , Blok-Tip L , Kaste D D . Designer drugs in herbal aphrodisiacs[J]. Forensic Science International, 2008, 177(2-3):0-0.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481