


Pyrithione sodium

Chemical Properties

CAS No. :	15922-78-8
Formula:	C ₅ H ₄ NNaOS
Molecular Weight:	149.15
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pyrithione sodium is isolated from cephalosporin molecule; some metal complexes of this have fumarate reductase inhibitory activity and may be useful against trypanosomes.
Targets(IC50)	Others,Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7047 mL	33.5233 mL	67.0466 mL
5 mM	1.3409 mL	6.7047 mL	13.4093 mL
10 mM	0.6705 mL	3.3523 mL	6.7047 mL
50 mM	0.1341 mL	0.6705 mL	1.3409 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Jing M, Chen Z, Li Z, Li F, Chen M, Zhou M, He B, Chen L, Hou Z, Chen X. Facile Synthesis of ZnS/N,S Co-doped Carbon Composite from Zinc Metal Complex for High-Performance Sodium-Ion Batteries. *ACS Appl Mater Interfaces*. 2018 Jan 10;10(1):704-712. doi: 10.1021/acsami.7b15659. Epub 2017 Dec 28. PubMed PMID: 29243910.
- Greenhalgh R, Greenhalgh M, Alshareef F, Robson GD. Application of green fluorescent protein to measure antimicrobial efficacy and the kinetics of cell death against *Escherichia coli*. *J Microbiol Methods*. 2017 Oct;141:67-72. doi: 10.1016/j.mimet.2017.08.005. Epub 2017 Aug 10. PubMed PMID: 28802720.
- Tuncay E, Turan B. Intracellular Zn(2+) Increase in Cardiomyocytes Induces both Electrical and Mechanical Dysfunction in Heart via Endogenous Generation of Reactive Nitrogen Species. *Biol Trace Elem Res*. 2016 Feb;169(2):294-302. doi: 10.1007/s12011-015-0423-3. Epub 2015 Jul 4. PubMed PMID: 26138011.
- Knaack JS, Zhou Y, Magnuson M, Silvestri E, Johnson RC. Performance of a novel high throughput method for the determination of VX in drinking water samples. *Anal Chem*. 2013 Mar 5;85(5):2611-6. doi: 10.1021/ac3036102. Epub 2013 Feb 12. PubMed PMID: 23402290.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481