


## Sachaconitine

## Chemical Properties

CAS No. :	1361-02-0
Formula:	C <sub>23</sub> H <sub>37</sub> N <sub>1</sub> O <sub>4</sub>
Molecular Weight:	391.552
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sachaconitine is a natural product for research related to life sciences. The catalog number is TN6341 and the CAS number is 1361-02-0.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.554 mL	12.7698 mL	25.5395 mL
5 mM	0.5108 mL	2.554 mL	5.1079 mL
10 mM	0.2554 mL	1.277 mL	2.554 mL
50 mM	0.0511 mL	0.2554 mL	0.5108 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481