


3',4'-Dihydroinfectopyrone

Chemical Properties

CAS No. : 1083198-57-5
 Formula: C₁₄H₁₈O₅
 Molecular Weight: 266.293
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	3',4'-Dihydroinfectopyrone is a useful organic compound for research related to life sciences. The catalog number is T126433 and the CAS number is 1083198-57-5.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7553 mL	18.7765 mL	37.553 mL
5 mM	0.7511 mL	3.7553 mL	7.5106 mL
10 mM	0.3755 mL	1.8777 mL	3.7553 mL
50 mM	0.0751 mL	0.3755 mL	0.7511 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481