


Chitin synthase inhibitor 3

Chemical Properties

CAS No. :	2416338-26-4
Formula:	C ₂₀ H ₁₉ N ₃ O ₄
Molecular Weight:	365.38
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Chitin Synthase Inhibitor 3 (Compound 2d) is a potent inhibitor of chitin synthase, exhibiting an IC ₅₀ value of 0.16 mM and demonstrating antifungal activity with a minimum inhibitory concentration (MIC) of 1 µg/mL against <i>Candida albicans</i> [1].
Targets(IC ₅₀)	Others,Antifungal

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7369 mL	13.6844 mL	27.3688 mL
5 mM	0.5474 mL	2.7369 mL	5.4738 mL
10 mM	0.2737 mL	1.3684 mL	2.7369 mL
50 mM	0.0547 mL	0.2737 mL	0.5474 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481