

BML-210

## Chemical Properties


CAS No. : 537034-17-6

Formula: C<sub>20</sub>H<sub>25</sub>N<sub>3</sub>O<sub>2</sub>

Molecular Weight: 339.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	BML-210 (CAY10433) is a new-type HDAC inhibitor.
Targets(IC50)	Apoptosis,HDAC

## Solubility Information

Solubility	DMSO: 60 mg/mL (176.77 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+90% Corn Oil: 2.5 mg/mL (7.37 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9461 mL	14.7306 mL	29.4612 mL
5 mM	0.5892 mL	2.9461 mL	5.8922 mL
10 mM	0.2946 mL	1.4731 mL	2.9461 mL
50 mM	0.0589 mL	0.2946 mL	0.5892 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Jayathilaka N, et al. Nucleic Acids Res. 2012 Jul;40(12):5378-88.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481