


Azido-PEG2-t-butyl ester

Chemical Properties

CAS No. :	20844-48-8
Formula:	C ₁₁ H ₂₁ N ₃ O ₄
Molecular Weight:	259.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Azido-PEG2-t-butyl ester is a PEG derivative containing a t-butyl ester and an azide group. The hydrophilic PEG spacer increases solubility in aqueous media. Under acidic conditions, the t-butyl protected carboxyl group can be deprotected.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8564 mL	19.2819 mL	38.5639 mL
5 mM	0.7713 mL	3.8564 mL	7.7128 mL
10 mM	0.3856 mL	1.9282 mL	3.8564 mL
50 mM	0.0771 mL	0.3856 mL	0.7713 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Sano K, Nakajima T, Miyazaki K, Ohuchi Y, Ikegami T, Choyke PL, Kobayashi H. Short PEG-linkers improve the performance of targeted, activatable monoclonal antibody-indocyanine green optical imaging probes. *Bioconjug Chem.* 2013 May 15;24(5):811-6. doi: 10.1021/bc400050k. Epub 2013 May 3. PubMed PMID: 23600922; PubMed Central PMCID: PMC3674550.

Harrison E, Coulter JA, Dixon D. Gold nanoparticle surface functionalization: mixed monolayer versus hetero bifunctional peg linker. *Nanomedicine (Lond).* 2016 Apr;11(7):851-65. Review. PubMed PMID: 27021417.

Augusto MT, Hollmann A, Porotto M, Moscona A, Santos NC. Antiviral Lipopeptide-Cell Membrane Interaction Is Influenced by PEG Linker Length. *Molecules.* 2017 Jul 15;22(7). pii: E1190. doi: 10.3390/molecules22071190.

PubMed PMID: 28714870; PubMed Central PMCID: PMC5776016.

Tuma R, Russell M, Rosendahl M, Thomas GJ Jr. Solution conformation of the extracellular domain of the human tumor necrosis factor receptor probed by Raman and UV-resonance Raman spectroscopy: structural effects of an engineered PEG linker. *Biochemistry.* 1995 Nov 21;34(46):15150-6. PubMed PMID: 7578129.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481