

Furanofukinin

Chemical Properties


CAS No. :

Formula: C₁₆H₂₄O₂

Molecular Weight: 248.366

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Furanofukinin is a useful organic compound for research related to life sciences and the catalog number is T125001.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0263 mL	20.1313 mL	40.2625 mL
5 mM	0.8053 mL	4.0263 mL	8.0525 mL
10 mM	0.4026 mL	2.0131 mL	4.0263 mL
50 mM	0.0805 mL	0.4026 mL	0.8053 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481