


Decaethylene glycol

Chemical Properties

CAS No. : 5579-66-8

Formula: C₂₀H₄₂O₁₁

Molecular Weight: 458.541



Keep away from moisture, Store at low temperature

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Decaethylene glycol (HO-PEG10-OH) is a PROTAC linker belonging to the PEG class and can be used to synthesize PROTAC molecules.
Targets(IC50)	PROTAC Linker

Solubility Information

Solubility	DMSO: 80.00 mg/mL (174.47 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 3.30 mg/mL (7.20 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1808 mL	10.9042 mL	21.8083 mL
5 mM	0.4362 mL	2.1808 mL	4.3617 mL
10 mM	0.2181 mL	1.0904 mL	2.1808 mL
50 mM	0.0436 mL	0.2181 mL	0.4362 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481