

Desirudin

Chemical Properties

CAS No. : 120993-53-5

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Desirudin

Biological Description

Description	Desirudin (CGP 39393) is a thrombin inhibitor that prevents thrombosis and venous stasis thrombosis. It is utilized for studying thrombocytopenia or platelet dysfunction associated with thrombocytopenia.
Targets(IC50)	Others,Thrombin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481