

TAN 1446C

## Chemical Properties


CAS No. : 148084-28-0

Formula: C17H24O5

Molecular Weight: 308.374

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	TAN 1446C is a natural product that can be used as a reference standard. The CAS number of TAN 1446C is 148084-28-0.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2429 mL	16.2143 mL	32.4286 mL
5 mM	0.6486 mL	3.2429 mL	6.4857 mL
10 mM	0.3243 mL	1.6214 mL	3.2429 mL
50 mM	0.0649 mL	0.3243 mL	0.6486 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481