

Actriol

Chemical Properties


CAS No. : 547-81-9

Formula: C₁₈H₂₄O₃

Molecular Weight: 288.38

Storage: Keep away from moisture, Store at low temperature
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Actriol (Epiestriol) is always used as Glucocorticoid Receptor agonist.
Targets(IC50)	Glucocorticoid Receptor, Antibacterial

Solubility Information

Solubility	DMSO: 55 mg/mL (190.72 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4676 mL	17.3382 mL	34.6765 mL
5 mM	0.6935 mL	3.4676 mL	6.9353 mL
10 mM	0.3468 mL	1.7338 mL	3.4676 mL
50 mM	0.0694 mL	0.3468 mL	0.6935 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Miller E, et al. J Interferon Cytokine Res. 1998 Nov;18(11):921-5.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481