

GSK-598809 L-tartrate

Chemical Properties


CAS No. : 1008521-08-1

Formula: C₂₆H₂₉F₄N₅O₇S

Molecular Weight: 631.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	GSK-598809 L-tartrate is an acute D3 antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5833 mL	7.9164 mL	15.8328 mL
5 mM	0.3167 mL	1.5833 mL	3.1666 mL
10 mM	0.1583 mL	0.7916 mL	1.5833 mL
50 mM	0.0317 mL	0.1583 mL	0.3167 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481