


F-050

Chemical Properties

CAS No. : 186835-06-3
 Formula: C₂₁H₂₀O
 Molecular Weight: 288.38
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	F-050 is a biochemical that has been shown to inhibit platelet aggregation induced by CaCl ₂ , arachidonic acid, collagen, adenosine diphosphate (ADP) and thrombin in guinea pigs, rabbits and rats in vitro. F-050 may be useful in the treatment of thrombotic diseases.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4676 mL	17.3382 mL	34.6765 mL
5 mM	0.6935 mL	3.4676 mL	6.9353 mL
10 mM	0.3468 mL	1.7338 mL	3.4676 mL
50 mM	0.0694 mL	0.3468 mL	0.6935 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481