


(Iso)-FK-480

## Chemical Properties

CAS No. : 167820-10-2  
 Formula: C<sub>26</sub>H<sub>19</sub>FN<sub>4</sub>O<sub>2</sub>  
 Molecular Weight: 438.45  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(Iso)-FK-480 is a novel orally available cholecystikinin A (CCK-A) antagonist for investigational treatment of chronic pancreatitis.
Targets(IC50)	Cholecystikinin Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2808 mL	11.4038 mL	22.8076 mL
5 mM	0.4562 mL	2.2808 mL	4.5615 mL
10 mM	0.2281 mL	1.1404 mL	2.2808 mL
50 mM	0.0456 mL	0.2281 mL	0.4562 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Satoh Y, et al. Studies on a novel, potent and orally effective cholecystikinin A antagonist, FK-480. Synthesis and structure-activity relationships of FK-480 and related compounds. Chem Pharm Bull (Tokyo). 1994 Oct;42(10): 2071-83.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481