

RR-src

Chemical Properties

CAS No. : 81156-93-6


Formula: C64H106N22O21

Molecular Weight: 1519.66

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tyrosine kinase substrate peptide
Targets(IC50)	Others,Src

Solubility Information

Solubility	20% acetonitrile / H2O: 1 mg/mL (0.66 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.658 mL	3.2902 mL	6.5804 mL
5 mM	0.1316 mL	0.658 mL	1.3161 mL
10 mM	0.0658 mL	0.329 mL	0.658 mL
50 mM	0.0132 mL	0.0658 mL	0.1316 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481