


H2N-PEG8-Hydrazide

Chemical Properties

CAS No. :

Formula: C19H41N3O9

Molecular Weight: 455.54



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	H2N-PEG8-Hydrazide is a PEG-based linker for PROTACs that joins two essential ligands, crucial for forming PROTAC molecules, enabling selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands joined by a linker: one targeting an E3 ubiquitin ligase and the other the target protein. By utilizing the intracellular ubiquitin-proteasome system, PROTACs selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1952 mL	10.976 mL	21.952 mL
5 mM	0.439 mL	2.1952 mL	4.3904 mL
10 mM	0.2195 mL	1.0976 mL	2.1952 mL
50 mM	0.0439 mL	0.2195 mL	0.439 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481