


3-(2-CHLORO-4-FLUORO-PHENYL)-PROPIONIC ACID

Chemical Properties

CAS No. :	174603-38-4
Formula:	C ₉ H ₈ ClFO ₂
Molecular Weight:	202.61
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-(2-CHLORO-4-FLUORO-PHENYL)-PROPIONIC ACID, with CAS No. 174603-38-4, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 3-(2-CHLORO-4-FLUORO-PHENYL)-PROPIONIC ACID provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9356 mL	24.678 mL	49.3559 mL
5 mM	0.9871 mL	4.9356 mL	9.8712 mL
10 mM	0.4936 mL	2.4678 mL	4.9356 mL
50 mM	0.0987 mL	0.4936 mL	0.9871 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481