


Boc-TIC-OH

Chemical Properties

CAS No. : 78879-20-6
 Formula: C₁₅H₁₉NO₄
 Molecular Weight: 277.32
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Boc-TIC-OH, with CAS No. 78879-20-6, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Boc-TIC-OH provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6059 mL	18.0297 mL	36.0594 mL
5 mM	0.7212 mL	3.6059 mL	7.2119 mL
10 mM	0.3606 mL	1.803 mL	3.6059 mL
50 mM	0.0721 mL	0.3606 mL	0.7212 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481