


Wnt pathway activator 1

Chemical Properties

CAS No. :	1360540-81-3
Formula:	C18H16O4
Molecular Weight:	296.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Wnt pathway activator 1 is a potent Wnt activator with an IC50 of 28-29 nM. Wnt pathway activator 1 is a potent Wnt activator for the treatment of androgenetic alopecia (AGA).
Targets(IC50)	Wnt/beta-catenin

Solubility Information

Solubility	DMSO: 60 mg/mL (202.48 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween-80+45% Saline: 2.5 mg/mL (8.44 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3747 mL	16.8737 mL	33.7473 mL
5 mM	0.6749 mL	3.3747 mL	6.7495 mL
10 mM	0.3375 mL	1.6874 mL	3.3747 mL
50 mM	0.0675 mL	0.3375 mL	0.6749 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Charlene F. Barroga, et al. Diketones and hydroxyketones as catenin signaling pathway activators. WO2012024404A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481