


2-(2-Chlorophenoxy)ethylamine

Chemical Properties

CAS No. :	26378-53-0
Formula:	C ₈ H ₁₀ ClNO
Molecular Weight:	171.624
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-(2-Chlorophenoxy)ethylamine, with CAS No. 26378-53-0, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2-(2-Chlorophenoxy)ethylamine provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8268 mL	29.1341 mL	58.2683 mL
5 mM	1.1654 mL	5.8268 mL	11.6537 mL
10 mM	0.5827 mL	2.9134 mL	5.8268 mL
50 mM	0.1165 mL	0.5827 mL	1.1654 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481