


5-(Azidomethyl)arauridine

Chemical Properties

CAS No. : 34349-03-6
 Formula: C₁₀H₁₃N₅O₆
 Molecular Weight: 299.24
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Nucleoside Derivatives - Azido-nucleosides, 5-Modified pyrimidine nucleosides, Arabino-nucleosides
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3418 mL	16.709 mL	33.418 mL
5 mM	0.6684 mL	3.3418 mL	6.6836 mL
10 mM	0.3342 mL	1.6709 mL	3.3418 mL
50 mM	0.0668 mL	0.3342 mL	0.6684 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481