

m-PEG6-Hydrazide

Chemical Properties

CAS No. : 1449390-64-0


Formula: C14H30N2O7

Molecular Weight: 338.4

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	m-PEG6-Hydrazide, a PEG-based linker for PROTACs, connects two essential ligands vital for the formation of PROTAC molecules, enabling selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands joined by a linker: one ligand targets an E3 ubiquitin ligase and the other targets the protein of interest. These compounds leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9551 mL	14.7754 mL	29.5508 mL
5 mM	0.591 mL	2.9551 mL	5.9102 mL
10 mM	0.2955 mL	1.4775 mL	2.9551 mL
50 mM	0.0591 mL	0.2955 mL	0.591 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481