


2'-Deoxy-N6-Phenoxyacetyladenosine

Chemical Properties

CAS No. :	110522-74-2
Formula:	C18H19N5O5
Molecular Weight:	385.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2'-Deoxynucleoside,2/6/8-Modified nucleoside
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5949 mL	12.9745 mL	25.9491 mL
5 mM	0.519 mL	2.5949 mL	5.1898 mL
10 mM	0.2595 mL	1.2975 mL	2.5949 mL
50 mM	0.0519 mL	0.2595 mL	0.519 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481