

## Compound T128287

## Chemical Properties


CAS No. :

Formula: C41H66O14

Molecular Weight: 782.97

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl (1R,4aS,6aS,6bR,9R,10S,12aR,14bS)-1-hydroxy-9-(hydroxymethyl)-1,2,6a,6b,9,12a-hexamethyl-10-(((2S,3R,4S,5S)-3,4,5-trihydroxyoxan-2-yl]oxy)-1,2,3,4,4a,5,6,6a,6b,7,8,8a,9,10,11,12,12a,12b,13,14-b-icosahydricene-4a-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2772 mL	6.3859 mL	12.7719 mL
5 mM	0.2554 mL	1.2772 mL	2.5544 mL
10 mM	0.1277 mL	0.6386 mL	1.2772 mL
50 mM	0.0255 mL	0.1277 mL	0.2554 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481