

(S)-2-Amino-3-cyclopropylpropanoic acid

Chemical Properties

CAS No. :	102735-53-5
Formula:	C ₆ H ₁₁ NO ₂
Molecular Weight:	129.157
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	(S)-2-Amino-3-cyclopropylpropanoic acid is an alanine derivative applicable for biochemical experiments and drug synthesis research.
Targets(IC50)	Amino Acids and Derivatives

Solubility Information

Solubility	H ₂ O: 80.00 mg/mL (619.40 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.7423 mL	38.7117 mL	77.4234 mL
5 mM	1.5485 mL	7.7423 mL	15.4847 mL
10 mM	0.7742 mL	3.8712 mL	7.7423 mL
50 mM	0.1548 mL	0.7742 mL	1.5485 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Luckose F, et al. Effects of amino acid derivatives on physical, mental, and physiological activities. Crit Rev Food Sci Nutr. 2015;55(13):1793-807.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481