


PPHP

Chemical Properties

CAS No. :	87864-20-8
Formula:	C ₁₁ H ₁₄ O ₂
Molecular Weight:	178.23
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	PPHP is a substrate utilized for measuring peroxidase enzymes and has been employed to quantify the peroxidase activity of COX-1 and COX-2.
Targets(IC50)	Others,COX

Solubility Information

Solubility	PBS (pH 7.2): >540 µg/mL,Sonication is recommended. DMSO: >7 mg/mL,Sonication is recommended. DMF: >10 mg/mL,Sonication is recommended. Ethanol: >10 mg/mL,Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6107 mL	28.0536 mL	56.1073 mL
5 mM	1.1221 mL	5.6107 mL	11.2215 mL
10 mM	0.5611 mL	2.8054 mL	5.6107 mL
50 mM	0.1122 mL	0.5611 mL	1.1221 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481