


2-(3-Chlorobenzamido)-3-(2-oxo-1,2-dihydroquinolin-4-yl)propanoic acid

Chemical Properties

CAS No. : 90098-05-8
 Formula: C₁₉H₁₅ClN₂O₄
 Molecular Weight: 370.79
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-(3-Chlorobenzamido)-3-(2-oxo-1,2-dihydroquinolin-4-yl)propanoic acid is a useful organic compound for research related to life sciences. The catalog number is T64461 and the CAS number is 90098-05-8.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6969 mL	13.4847 mL	26.9694 mL
5 mM	0.5394 mL	2.6969 mL	5.3939 mL
10 mM	0.2697 mL	1.3485 mL	2.6969 mL
50 mM	0.0539 mL	0.2697 mL	0.5394 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481