


Cycloartanol

Chemical Properties

CAS No. :	4657-58-3
Formula:	C ₃₀ H ₅₂ O
Molecular Weight:	428.745
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cycloartanol is a natural product from Citrus medica.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3324 mL	11.6618 mL	23.3236 mL
5 mM	0.4665 mL	2.3324 mL	4.6647 mL
10 mM	0.2332 mL	1.1662 mL	2.3324 mL
50 mM	0.0466 mL	0.2332 mL	0.4665 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481