


Selinidin

Chemical Properties

CAS No. :	19427-82-8
Formula:	C ₁₉ H ₂₀ O ₅
Molecular Weight:	328.364
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Selinidin decreases phosphorylation of phospholipase C-gamma1, p38 mitogen-activated protein kinase, and IkappaB-alpha upon FcepsilonRI stimulation.
Targets(IC50)	IkB/IKK,p38 MAPK,TNF

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0454 mL	15.2272 mL	30.4544 mL
5 mM	0.6091 mL	3.0454 mL	6.0909 mL
10 mM	0.3045 mL	1.5227 mL	3.0454 mL
50 mM	0.0609 mL	0.3045 mL	0.6091 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481