


Delphinidin-3-O-arabinoside chloride

Chemical Properties

CAS No. : 171370-55-1
 Formula: C₂₀H₁₉ClO₁₁
 Molecular Weight: 470.81
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Delphinidin-3-O-arabinoside chloride is a dietary antioxidant.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.124 mL	10.620 mL	21.240 mL
5 mM	0.4248 mL	2.124 mL	4.248 mL
10 mM	0.2124 mL	1.062 mL	2.124 mL
50 mM	0.0425 mL	0.2124 mL	0.4248 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481