


## Segetoside K

## Chemical Properties

CAS No. :	261160-05-8
Formula:	C <sub>54</sub> H <sub>86</sub> O <sub>26</sub>
Molecular Weight:	1151.26
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Segetoside K is a useful organic compound for research related to life sciences. The catalog number is T125175 and the CAS number is 261160-05-8.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8686 mL	4.3431 mL	8.6861 mL
5 mM	0.1737 mL	0.8686 mL	1.7372 mL
10 mM	0.0869 mL	0.4343 mL	0.8686 mL
50 mM	0.0174 mL	0.0869 mL	0.1737 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481