


L-156,602

Chemical Properties

CAS No. : 125228-51-5
 Formula: C₃₈H₆₄N₈O₁₃
 Molecular Weight: 840.973
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	L-156,602 has a wide range of applications in life science related research.
Targets(IC50)	Others, Complement System

Solubility Information

Solubility	Ethanol: Soluble DMSO: Soluble DMF: Soluble (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1891 mL	5.9455 mL	11.891 mL
5 mM	0.2378 mL	1.1891 mL	2.3782 mL
10 mM	0.1189 mL	0.5946 mL	1.1891 mL
50 mM	0.0238 mL	0.1189 mL	0.2378 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481